

RELIABLE ANALYSIS	INC.

Tel: (248)-588-9770 • Fax: (248)-588-9769 E-mail: klang@ralab.com

REPORT NUMBER	1101397R-16676
PURCHASE ORDER	11000922
TEST DATE	8/5/11-8/22/11
REPORT DATE	8/22/11
TOTAL PAGE (S) NUMBER	1/6

REPORT FOR

Palziv 3745 Napier Rd. Plymouth, MI 48170

Attn: Mark Williams

Work Requested

Perform the following tests on materials submitted for testing in accordance with the instructions of the submitting engineer and the indicated laboratory procedures:

- 1. Fogging (GMW3235-B, B-Y1-T1)
- 2. Odor (GMW3205, Wet and Dry, 24 H @ 70 °C)
- 3. TVOC (GMW8081)
- 4. Mildew (GMW3259, 14d @ 40°C)

WSS-M99P32-B

3.11 Fogging (SAE J1756, 100°C heat, 21°C cool, 1h and 16h post)

Sample Description

Four materials, submitted for testing on August 5, 2011, were identified as:

- 1. Vizion 2 PCF-V 2.0
- 2. Vizion 3 PCF-V 3.0
- 3. Vizion 4 PCF-V 4.0
- 4. Vizion 2 PCF Soft- V 2.0S, Grey w/ EVA

Work Performed

The specimens were conditioned for a minimum of 24 hours at the standard laboratory conditions of 23 ± 2 °C and $50 \pm 5\%$ R.H. prior to testing. The tests were performed in accordance with the above indicated procedures.

Reliable Analysis Inc. letters and reports are for the exclusive use of our customers to whom they are addressed and shall not be reproduced except in full without written approval. The letters and reports apply only to those samples tested and are not necessarily indicative of the qualities of apparently identical or similar products. Samples will be retained for a maximum of thirty days after the date of the final report. The name Reliable Analysis Inc. or its insignia or reports are not to be used by the customer in any forms of media advertising.



PAGE NUMBER

2/6

<u>Test Results</u>

1. Fogging (GMW3235B, Code 3, 18h @ 90°C, 21°C cool, 24h post)

The specimens were preconditioned in a desiccator filled with silica gel for 1 day. After conditioning, each specimen was put in a beaker, and the beaker was covered with a disk of aluminum foil and a glass plate. The beakers were then heated in a Hart Fog Test Chamber for 18 hours at 90 °C heating unit temperature and 21 °C cooling plate temperature. The aluminum foil disk was weighed before the test and 24 hours later after the test. Then the weight of condensate was determined (W_f - W_0). **Requirement:** NA

Sample ID	<u>W₀ (mg)</u>	<u>W_f (mg)</u>	Weight of Condensate (mg)
V 2.0			
Specimen 1	307.5	309.4	1.9
Specimen 2	307.3	309.2	1.9
Specimen 3	305.4	307.5	2.1
		Mean:	2.0
V 3.0			
Specimen 1	304.9	306.8	1.9
Specimen 2	308.5	310.0	1.5
Specimen 3	306.1	307.5	1.4
		Mean:	1.6
V 4.0			
Specimen 1	306.9	308.9	2.0
Specimen 2	309.2	311.0	1.8
Specimen 3	308.3	310.4	2.1
		Mean:	2.0

Odor (GMW3205, Wet and Dry, 24 H @ 70 °C) The specimens were tested at 70 °C for 24 hours. A ratio of specimen size to container size of 1:20 was used. Requirement: NA

			Par	nelist			Average
<u>Sample ID</u>	<u>A</u>	<u>B</u>	<u>C</u>	D	E	F	Rating
<u>Coding A- Wet</u> V 2.0							
Specimen 1	4	8	9	7	8	7	7.2
Specimen 2	4	8	9	7	8	7	7.2
Specimen 3	4	8	9	7	8	7	7.2
Rating Scale:							
U			10	Odo	rless		
			8	Trac	e		
			6	Tole	rable		
			4	Ann	oying		
			2	Seve			



PAGE NUMBER

3/6

Test Results (continued)

 Odor (GMW3205, Wet and Dry, 24 H @ 70 °C) The specimens were tested at 70 °C for 24 hours. A ratio of specimen size to container size of 1:20 was used.
 Requirement: NA

Panelist Average E F Rating Sample ID A B С D Coding A-Wet V 3.0 Specimen 1 7.7 Specimen 2 7.7 Specimen 3 7.7 V 4.0 Specimen 1 8.8 Specimen 2 8.8 Specimen 3 8.8 **Coding B- Dry** V 2.0 Specimen 1 Specimen 2 Specimen 3 V 3.0 Specimen 1 5.8 Specimen 2 5.8 Specimen 3 5.8 V 4.0 Specimen 1 Specimen 2 Specimen 3

3. TVOC (GMW 8081)

The test conditions were as listed below:

Parameters of headspace feeding flask: Temperature: 120°C oven, 125°C loop and 135°C transfer line Time: GC cycle time 60 min
Equilibration Time: 60 min
Pressurize Time: 0.1min
Loop Fill Time: 0.1 min
Loop Equilibration Time: 0.0 min
Inject Time: 1.0 min



PAGE NUMBER

Test Results (continued)

3. TVOC (GMW 8081)

Measuring condition:

GC furnace temperature: 4 min isothermal condition at 50°C, heating to 300°C at a rate of 7°C/min and then 300°C for 5 min;

Pressure before entering into the column: 10 psi.

The "Relative Emission" was calculated in accordance with the following formula:

Relative Emission = (Total Area- Air Area) / Air area

Requirement: NA

Sample ID	<u>A</u> sum	<u>A_{air}</u>	A Emission	Relative Emission
V 2.0				
Specimen 1	380665182	371261777	9403405	0.0253
Specimen 2	467301033	451566935	15734098	0.0348
Mean				0.0301
V 3.0				
Specimen 1	342027922	333286985	8740937	0.0262
Specimen 2	348902261	338479465	10422796	0.0308
Specimen 3	383487056	372801917	10685139	0.0287
Mean				0.0286
V 4.0				
Specimen 1	328153232	328153232	0	0.000
Specimen 2	320580356	320580356	0	0.000
Specimen 3	321153972	3211539972	0	0.000
Mean				0.000

4. Mildew (GMW3259, 14d @ 40°C)

The 2"X 2" specimens were each placed in a sealed glass container with enough water to last the entire 14 day test period. The specimens were suspended above the water and not in contact with the glass. After 14 days in the incubator @ 40° C the containers were removed and allowed to cool. The odor was evaluated first and then each specimen was visually evaluated for the presence of mildew. **Requirement:** NA

Sample ID	<u>Odor</u>	Visual Examination	
V 2.0			
Specimen 1	No offensive odor	No noticeable growth	
Specimen 2	No offensive odor	No noticeable growth	
Specimen 3	No offensive odor	No noticeable growth	
V 3.0			
Specimen 1	No offensive odor	No noticeable growth	
Specimen 2	No offensive odor	No noticeable growth	
Specimen 3	No offensive odor	No noticeable growth	



REPORT NUMBER

PAGE NUMBER

Test Results (continued)

4. Mildew (GMW3259, 14d @ 40°C)

<u>Sample ID</u>	<u>Odor</u>	<u>Visual Examination</u>
V 4.0		
Specimen 1	No offensive odor	No noticeable growth
Specimen 2	No offensive odor	No noticeable growth
Specimen 3	No offensive odor	No noticeable growth

WSS-M99P32-B

3.11 Fogging (SAE J1756, 100°C heat, 21°C cool, 1h and 16h post **Requirement:** Minimum fog #70. No formation of clear film, droplets or crystals

Sample ID	<u>R</u> 0	<u>R</u> f	Fog [#]	Visual Examination	Evaluation
	<u> </u>				
V 2.0S					
Specimen 1	151.0	-	-	Large, wet droplets. Cannot read	Did not meet requirement
Specimen 2	151.0	-	-	Large, wet droplets. Cannot read	Did not meet requirement
Specimen 3	151.0	-	-	Large, wet droplets. Cannot read	Did not meet requirement
Mean:					
V 2.0S					
Specimen 1	151.0	115.0	76.2	Droplets present	Did not meet requirement
Specimen 2	151.0	117.0	77.5	Droplets present	Did not meet requirement
Specimen 3	151.0	117.0	77.5	Droplets present	Did not meet requirement
Mean:			77.0		

Test Equipment

- 1. E-070, Precision Scientific Oven (Calibration Due: 1/12)
- 2. E-103, Denver Analytical Balance (Calibration Due: 4/12)
- 3. E-005, Hart Scientific Fog-Test Chamber (S/N: 125, Calibration Due: 7/12)
- 4. E-055, Incubator (Calibration Due: 1/12)
- 5. E-006, Hart Scientific Fog-Test Chamber (S/N: 140, Calibration Due: 7/12)
- 6. E-007, Glossmeter (Calibration Due: 8/11)



REPORT NUMBER

PAGE NUMBER

6/6

Sample Disposition

The tested samples are being returned with this report for your further evaluation.

Reliable Analysis, Inc.

Fen - Jen Day

Ken-jen Lang President

SH/ethompson/reports2011